

COMPO-JECT™

Easy to use. Ergonomic. Built to last.

COMPO-JECT™ is a lightweight, ergonomic composite dispensing gun designed for smooth, controlled handling in everyday restorative procedures. Its balanced design and proven functionality provide comfortable and effortless operation.



Made from high grade materials.

- **Takes universal** pre-filled unit-dose tips.
- **The tip** can be rotated in all directions – and is safely fixed.



Easy and quick replacement of the tip.

Effortless dispensing. Gentle on hand and wrist.

	Weight	Width	Height	Depth
REF 20006	111 g	125 mm	97,5 mm	17,3 mm

Materials: High grade stainless steel and glass-reinforced polymer.

Withstands therm disinfection and moist heat sterilization (autoclave).

Made in Denmark



Scan and watch video

Version 01



Manufactured by
RØNVIG
Dental Mfg. A/S
www.ronvig.com
contact@ronvig.com